# STARFIT™

## ADULT NASOGASTRIC FEEDING TUBE

Single-use device.

Length of Use up to 30 days.

Sterile device(E0).

Do not reuse, do not resterilize.

Does not contain Latex or DEHP or BPA.

Store dry and protected from light.

Manufacturer :

**HOIST MEDICAL**

23 Quai de Rive Neuve

13007 Marseille France

Tel : + 33 (0) 4 91 76 73 29

hoist-medical.com

**EN**
**FR**
**DE**
**IT**
**ES**

### Indications

STARFIT™ tubes are designed for nasogastric enteral feeding, passive gastric drainage/aspiration or nasointestinal feeding in patients requiring enteral nutrition or NG/NI access.

### Contraindications

Organic intestinal occlusion.

Repeated or incoercible vomiting.

High flow fistulas.

Severe acute anchoratitis in the first 8 days.

Severe diarrheas.

Insufficient effective intestinal absorbent surface.

### Description


The tubes are made of medical grade polyurethane with a 40% barium radiopaque line.

The NGG, NGW, NGGV, NGGBY and NGWBY references are equipped with a stainless steel stylet that facilitates insertion.

Weighted NGW and NGWBY tubes have a distal end weighted with tungsten.

STARFIT™ tubes are available in size Fr 06 Fr 08; Fr 10; Fr 12 and Fr 14 (Fr 16 and Fr 18 on request). The length is 80 cm 92 cm 120cm 140 cm.

The numeric centimetric markings printed on the tube facilitate insertion, documentation and allow for migration monitoring. All feeding tubes have ENFit® connector in accordance with ISO8036 9-3.

### Precautions and warnings

To be used by trained medical staff, trained patients or trained patient family.

Do not use if package is opened or damaged.

Check the expiry date and do not use after this date.

If lubrication is deemed necessary, use water based lubricant, do not use mineral oil or petroleum jelly.

Always check the position of the tube:

By aspiration of gastric fluid (with a 50 or 60 mL ENFit® enteral feeding syringe).

The position of the feeding tube on high-risk patients (comatose patients, non-reactive patients, etc.) may be confirmed by radiology.

### Insertion of the tube into the trachea, bronchi or lungs creates serious risks for the patient

Check that the tube has not wrapped into the patient's back of the throat. Never reintroduce the stylet when the tube is already in place.

Check the tube position periodically and systematically before each administration. Document the cm marking at patients' nares to observe for signs of migration.

The use of a syringe with a volume of less than 50 mL for positioning verification and flushing should be used with care as may generate pressures greater than the tube rupture pressure. Small volume syringes 1 mL to 20 mL are advised for drug administration only, with gentle pressure.

The tube should be flushed with 20 to 30 mL of water (in a 50 or 60 mL ENFit® syringe) whenever feeding is interrupted and systematically every 6 hours.

In case of obstruction

Do not attempt to unblock the tube, replacement is strongly recommended, especially if tube position cannot be confirmed. The use of a stylet to unblock the tube is prohibited.


The device to be destroyed after use according to the facility's biological waste management protocol. For home use, health workers will take care of waste management. Do not dispose of the device in household waste.

### Instructions for use

Please note: The following protocol is given for information purposes only, special situations may require certain modifications. Refer to local policies

Explain the following technique to the patient:

- Wash your hands before handling the device and its components and follow the hygiene rules established by the facility.
  - Place the patient in a seated position, slightly tilted back.
  - Determine the length of the tube required for it to enter the stomach, the distance from the earlobe to the nostril and the distance from the nostril to the xiphoid appendix. Use the centimeter graduations printed on the tube from the proximal end to avoid introducing it over an excessive length.
  - The tubes are pre-lubricated and do not require any special preparation before insertion. Lubrication can be improved by coating the tube with a water-based lubricant.
  - Insert the tube through one of the nostrils. Advance the tube through the esophagus into the stomach. Swallowing movements of the patient with or without fluid absorption facilitate the passage of the tube.
- Note: if the patient starts coughing during this phase, the tube has probably been inserted into the trachea, remove the tube immediately and repeat the procedure. In case of resistance, do not insist, remove the tube completely and repeat the procedure.
- For a tube equipped with a stylet:
- Before removing the stylet, always confirm the correct position of the tube (see precaution for use).
  - Remove the stylet and discard.
  - Tape the tube to the patient's nose and/or cheek if the chosen feeding site is the stomach.
  - For feeding in the duodenum or jejunum, leave the necessary length free (About 20 to 25 cm) before attaching the tube. (image below)

The tip of the tube will be driven by intestinal peristalsis. STARFIT™ tungsten weighted tubes are more suitable in this case. Check the progress of the tube every hour by checking the depth marks. Confirm tube position by local, national protocols/guidelines prior to feeding, administration of medications or if you have any concerns regarding tube position.

### Medication Administration

Medication administration should be separated from nutrition.

Any interaction between nutrients and medication must be approved by a physician or pharmacist. Check the prescription to ensure that the medication should be on an empty or full stomach. Administer medication in liquid form, to avoid clogging. Consult a physician or pharmacist to determine if it is possible to dilute it before use. If certain medications are only available in tablet form, consult a physician or pharmacist before crushing and mixing the tablet with water.

The consecutive administration of several medications should be undertaken on the advice of a physician or pharmacist with flushing between each medication.

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