STARFIT™

NASOENTERIC FEEDING TUBE

Single-use device

Length of Use up to 30 days

Sterile device (EO).

Do not reuse and do not resterilize.

Does not contain Latex or DEHP or BPA.

Store dry and protected from light.

Manufacturer :

HOIST MEDICAL

23 Quai de Rive Neuve

13007 Marseille France

Tel : + 33 (0) 4 91 76 73 29

hoist-medical.com

0297

EN

FR

DE

IT

ES

SE

Indications

STARFIT™ nasogastric enteral tube designed for continuous or intermittent enteral feeding and gastric access of neonate and paediatric patients.

Contraindications

Mechanical occlusion of the small bowel or colon (ileus).

Repeated or incoercible vomiting.

Active gastrointestinal hemorrhage.


Problems of coagulation.

Description

Clear medical grade polyurethane feeding tube with a 40% barium radiopaque line.

Numeric centimeter markings printed on the tube to facilitate insertion, documentation and allow for migration monitoring.

All feeding tubes have ENFit® connector in accordance with ISO80369-3 standard.

Precautions and warnings

To be used by trained medical staff, trained patient or trained patient family.

Do not use if package has been opened or damaged.

Check the expiry date and do not use after this date.

If lubrication is deemed necessary, use water based lubricant do not use mineral oil or petroleum jelly.

Always check the position of the tube:

By aspiration of gastric fluid (with a 50 or 60 mL ENFit® enteral syringe)

The position of the feeding tube on high-risk patients may be confirmed by radiology.

Insertion of the tube into the trachea, bronchi or lungs creates serious risks for the patient

Check that the tube has not wrapped into the patient's back of the throat.

Check the tube position periodically and systematically before each administration. Document the cm marking at patients' nares to observe for signs of migration.

The use of a syringe of less than 50 mL for positioning checking and flushing may result in pressures higher than the tube rupture pressure. The use of a syringe of small volume 1 ml to 20 ml are advised for the administration of drugs, with gentle pressure.

The tube should be flushed with a volume of water equal to approximately 3 to 4 times the internal volume of the tube, see table below (with a 50 or 60 mL ENFit® syringe), whenever feeding is interrupted and systematically every 6 hours.

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Interior volume of STARFIT™ feeding tubes.

In case of obstruction, do not attempt to unblock the tube, replacement is strongly recommended, especially if the tube tip position cannot be confirmed. The use of a stylet to unblock the tube is prohibited.

The device to be destroyed after use according to the facility's biological waste management protocol. For home use, health care workers will take care of waste management. Do not dispose of the device in household waste.

Instructions for use

Please note: the following protocol is given for information purposes only, special situations may necessitate certain modifications. Refer to local policies.

- Wash your hands before handling the device and its components and follow the hygiene rules established by the facility.

- Determine the length of the tube required for it to enter the stomach, the distance from the earlobe to the nostril and the distance from the nostril to the xiphoid appendix. Use the centimeter graduations printed on the tube to avoid introducing it over an excessive length.

- Place the patient in a seated position, slightly tilted back.

- Insert the tube through one of the nostrils. Insert the tube through the esophagus into the stomach. Swallowing movements of the patient with or without fluid absorption facilitate the passage of the tube.

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