



Hoist
medical

Improving quality of life with assisted feeding

INTRODUCING THE FUTURE OF ENTERAL FEEDING

PURFIT™

THE FIRST AND ONLY POLYURETHANE BALLOON GASTROSTOMY BUTTON

"INNOVATING WITH EXCELLENCE, CARING WITH CONFIDENCE."

"Elevate Patient Comfort and Confidence with Our Advanced Polyurethane Balloon Technology - Designed for Reliability and Ease of Use."

Designed & Manufactured with Excellence

With over **35 years of experience**, specialising in Enteral Nutrition Devices, managing the complete process from concept, design to manufacturing and supply in our **ISO-certified facilities**. Our products comply and are certified to the latest **EUMDR standard**. Our highly experienced and committed team ensures **quality, safety, and innovation** in every product we create.

The First & Only Polyurethane Balloon Gastrostomy Low Profile Button

A Revolutionary Choice for Clinicians & Patients. Designed for patients requiring short and long term enteral nutrition from **Hospital to Home**.



PURFIT™
GASTROSTOMY BUTTON

Why Choose Polyurethane?

Unlike silicone balloons, our polyurethane balloon technology delivers superior reliability and ease of use, reducing the need for interventions while improving patient confidence and comfort.

- ✓ **Increased Durability & Longevity**
Polyurethane material offers enhanced wear resistance*.
- ✓ **Consistent Balloon Shape**
Provides confidence - no need to massage and reshape before insertion.
- ✓ **Less Interventions**
No need for weekly balloon checks.
- ✓ **Atraumatic Inverted Shaft Tip**
Reduces risk of gastric abrasion or trauma.
- ✓ **Accessible Angled Inflation Port**
Angled away from the abdomen for patient comfort and ease of use.
- ✓ **Advanced CLICK AND CONNECT™ Technology**
Secure and intuitive connection for extension sets, simple for limited dexterity.
- 📏 **Wide range of sizes from 12 Fr to 24 Fr**



Polyurethane PEG catheters seem to be more resistant to deterioration than silicone PEG catheters, and at present they should be preferred for long-term enteral feeding via PEG*.



Hoist medical

Improving quality of life with assisted feeding

Give your patients and carers
the confidence they deserve
with a gastrostomy button that
is designed for **ease of use,**
reliability, and superior
patient care.

EXPERIENCE THE DIFFERENCE!

**CONTACT US TODAY
TO LEARN MORE**



[hoist-medical.com](https://www.hoist-medical.com)

**PURFIT™
GASTROSTOMY BUTTON**



Sales and General Enquiries
+33(0) 4 91 76 73 29
commercial@hoist-medical.com