



**Hoist**  
medical

## NEONATE PRODUCT CATALOGUE

Nurturing the Future: Where Every Drop Matters

Advanced Enteral  
Feeding Solutions  
for Neonatal Care

[hoist-medical.com](http://hoist-medical.com)

## Introduction

---

Welcome to our comprehensive catalogue of enteral feeding devices designed specifically for neonates. In the delicate world of neonatal care, proper nutrition is not just a crucial factor; it is the cornerstone of survival, growth, and healthy development. Ensuring that our tiniest babies receive the essential nutrients they need can be challenging, yet it is of utmost importance to survive and thrive.

**Hoist Medicals'** collection of enteral feeding solutions is tailored to meet the unique nutritional requirements of neonates. We understand that every neonatal baby has different needs, and our range of products reflects this diversity, offering safe, reliable, and innovative solutions to support healthcare professionals in delivering the best possible care.

From the latest advancements in feeding tubes to pumps and accessories, each product in our catalogue is designed with the highest standards of safety and efficacy in mind. We are committed to supporting neonatal healthcare providers with solutions that enhance feeding efficiency, minimise risks, and promote the overall well-being of neonates, because we know Every Drop Matters.

Explore our catalogue to discover a range of devices that can help ensure your patients receive the precise nutrition and medication delivery they need.



# CONTENTS

---

STARFIT™ NASOGASTRIC FEEDING TUBES.....	5
STARFIT™ ENTERAL SYRINGE PUMP.....	6
STARFIT™ ENTERAL ACCESSORIES.....	7 - 11
PURFIT™ GASTROSTOMY FEEDING.....	12- 15
PRODUCT CODES.....	16 - 17
ORDERING INFORMATION.....	18
SUSTAINABILITY.....	19
CONTACT INFORMATION.....	20

## THE ROLE OF ENTERAL DEVICES IN NUTRITION DELIVERY

Enteral devices, including their connectors, dead space and priming volumes, all create opportunities to reduce volume losses of Human Breast Milk in the enteral delivery system to ensure that babies receive their optimal nutrition, including lipids and fats.

### LIPID DELIVERY CHALLENGES

Vertical orientation is, contraindicated in many enteral syringe pumps, the horizontal orientation and position of the syringe tip to 'off set' ensures a **lipids first** delivery. The same applies to bolus delivery.



### SAFETY & MISCONNECTION PREVENTION

The initiative to prevent misconnections between different fluid delivery systems within the healthcare system resulted in the creation of the ISO 80369 series of standards. The ISO 80369-3 standard describes the dimensions and functional requirements for enteral small-bore connectors known as ENFit®.

This universal connector also simplifies and standardises equipment to allow for a safer, NICU step down, transfers between wards and hospitals – ensuring the transfer of babies, avoids unnecessary change of feeding tubes and accessories.




**TOGETHER, LET'S MAKE A DIFFERENCE IN NEONATAL CARE**

## POLYURETHANE

### WITH 40% BARIUM RADIOPAQUE STRIPE

Polyurethane material provides a larger internal lumen than that of PVC or silicone, maximising flow and reducing risk of blockages. Polyurethane is preferred by clinicians for babies as it reduces kinking during insertion and softens at body temperature, ensuring comfort and better tolerance when in place.



- ✓ Printed NUMERIC cm markings to facilitate insertion, provide visual cue for potential tube migration and support good documentation
- ✓ Rounded tip for atraumatic insertion, maximising babies comfort
- ✓ Closed distal tip with side exit ports reducing risk of adherence to mucosal of the stomach, no dead space helps avoid risk of bacterial contamination
- ✓ Clear tubing with blue radiopaque stripe for clear visualisation on x-ray
- ✓ Single port with simple and secure press in tethered cap


Purple ENFit® connector identifies for Enteral Use Only

Extensive size range offering from 4 Fr, 50 cm to 10 Fr, 120 cm

Approximate nasogastric tube priming volumes			
Code	50 cm	80 cm	120 cm
NG04-XX	0,5 ml	0,8 ml	1,1 ml
NG05-XX	0,6 ml	1,0 ml	1,4 ml
NG06-XX	0,7 ml	1,1 ml	1,6 ml
NG08-XX	0,9 ml	1,4 ml	2,1 ml
NG10-XX	1,7 ml	2,8 ml	4,2 ml

## WE KNOW EVERY DROP MATTERS

The STARFIT™ Enteral Syringe Pump is engineered for efficient and precise nutrition delivery. It boasts a large, easy-to-read screen, an intuitive keyboard, an alarm beacon, and a familiar user interface for smooth operation. Specifically calibrated for STARFIT™ Enteral Syringes, it features distinctive purple labelling, ensuring quick and safe identification for enteral use only.




## EXTENSION SETS

Guidelines recommend extension sets are changed every 4 hours to reduce risk of bacterial growth and contamination<sup>1</sup>

ENFit® Extension Sets to connect an ENFit® device/pump to a feeding tube.

- ✓ PVC with purple stripe and connectors - easily identified for enteral
- ✓ Single use
- ✓ Sterile
- ✓ Available in 75 cm and 150 cm lengths
- ✓ Pinch clamp
- ✓ Tethered cap
- ✓ Individually packaged
- ✓ DEHP/Latex Free



<sup>1</sup> Guidelines for the Preparation and Handling of Expressed and Donor Breast Milk and Specialist Feeds for Infants and Children in Neonatal and Paediatric Health Care Settings

HOIST MEDICAL carefully design and manufacture a wide range of accessories to help with the delivery of Human Breast Milk (HBM), Enteral feeds and Medicines. We understand the crucial role that the preparation and administration of nutrients, HBM, and medications play in neonatal health, where Every Drop Matters, we are committed to supporting their growth through our specialised solutions.

### MILK STRAWS

For the drawing up of larger volumes or thicker preparations



Code	Size/Description	Unit of Supply
<b>MS01</b>	Milk Straw	50 per box
<b>MS01WB</b>	Milk Straw Wide Bore	50 per box

### ENFit® CAPS

For ENFit® connectors to allow for closed system when not in use, and to provide a closed system when preparing HBM, Medications in advance of delivery.



Code	Size/Description	Unit of Supply
<b>MC01</b>	Male Cap	100 per box
<b>MC50</b>	Male Cap x 50	1000 per box
<b>MFCAP</b>	Male Female Cap	100 per box



Male Cap Available by 1 or 50 single package.

Code: MC01, MC50

Male-Female Connector

Code: MFCAP

### ORAL MEDICINE DELIVERY - DOSE AID AND DOSE STRAW

Provides a comfortable solution for babies, when administrating medications orally to a Neonate or Paediatric patient, alleviating concerns over placing threaded syringe tips or Low Dose tabs directly into babies mouths.



**MEDICINE DOSE AID**

Code: MF01

50 per box




**MEDICINE DOSE STRAW**

Code: MFS01

100 per box

### EXTENSION SETS

To allow connection between feeding tube and enteral syringe pump or gravity syringe



Code	Size/Description	Unit of Supply
<b>NGE-75</b>	75 cm	40 per box
<b>NGE-150</b>	150 cm	25 per box

### SAMPLING STRAW



Code	Size/Description	Unit of Supply
<b>SS01</b>	Sampling Straw	100 per box

Additional accessories to support an aseptic technique to draw up medications and milk.

### BOTTLE ADAPTORS

Bottle adaptors for use with ENFit® syringes , to allow drawing up of medication straight to syringe, reducing risk of contamination and secure closure between every access.



Code	Size/Description	Unit of Supply
BA-L01	20-21,5 mm	100 per box
BA-M01	18,5-20,5 mm	100 per box
BA-S01	18-19,5 mm	100 per box
BA-XL01	25,5-27 mm	100 per box
BA-XS01	15-16,5 mm	100 per box
UBA01	Universal 14,5-27,5 mm	30 per box

### NEEDLE STRAW, STERILE

For use with multi-dose bottle with a rubber bung



Code	Size/Description	Unit of Supply
NS01	Needle Straw	100 per box

### RECTAL STRAW

For the delivery of rectal medication



Code	Size/Description	Unit of Supply
RK75	Rectal Straw	100 per box
RK150	Rectal Straw	50 per box

### MALE MALE TRANSFER CONNECTOR

ENFit Male/Male connector for medicines transfer between enteral ENFit® syringes



Code	Size/Description	Unit of Supply
<b>TR01</b>	Male Male Transfer Connector	100 per box

### 3 WAY STOPCOCK ENFit®

To allow a closed system to be maintained with feeding tube, during delivery of nutrition and medication delivery



Code	Size/Description	Unit of Supply
<b>SK01</b>	3 Way Stopcock	25 per box

### ENTERAL BRUSH



Code	Size/Description	Unit of Supply
<b>BR01</b>	Enteral Cleaning Brush	10 x 10 bags per box

### DRAINAGE BAG, EXTENSION TUBING AND HANGER

Allows a closed system for the collection and measurement of gastric contents using an ENFit® feeding device, vented to help relieve gastric pressure and prevent air build-up in bag



Code	Size/Description	Unit of Supply
<b>DBS-100</b>	100 ml	30 per box
<b>DBS-250</b>	250 ml	30 per box

## ENFit® ENTERAL SYRINGES

### MAXIMUM NUTRITION DELIVERY

An off-set tip is designed to deliver colostrum and lipids first when placed in a horizontal syringe pump, making the most of the available nutrition. Validated with STARFIT™ Enteral Syringe Pump.

The ENFit® connector offers maximum safety for the patient by avoiding connection/ wrong route errors.



- ✓ Clear Barrel with purple plunger, easy visual identification for enteral
- ✓ Clear printed measurement graduations for dose accuracy
- ✓ Low-Dose Tip on 1 ml, 2.5 ml & 5 ml syringes avoids an over or underdose
- ✓ Offset tip promotes delivery of colostrum and lipids - 20 ml & 60 ml
- ✓ "Enteral use" printed on the syringe barrel

STARFIT™ ENTERAL ENFit® SYRINGES - SINGLE USE		
Code	Size/Description	Unit of Supply
SE01LD	1 ml low dose	100 per box
SE025LD	2.5 ml low dose	100 per box
SE05LD	5 ml low dose	100 per box
SE10	10 ml	100 per box
SE20	20 ml	80 per box
SE60	60 ml	60 per box

## GASTROSTOMY

---

The choice of size and type of gastrostomy feeding tube will depend on a number of factors including the expected duration of feeding, the reason nutrition or access to the digestive tract is required, and the condition and anatomy of your patient.



Hoist Medical have innovated the simple silicone button and have developed the **only Polyurethane Balloon button**, providing the all the benefits of polyurethane PEG in a convenient and easy to use low profile button.

**'Polyurethane** PEG catheters seem to be more resistant to deterioration than silicone PEG catheters, **and at present they should be preferred for long-term enteral feeding** via PEG'.<sup>1</sup>

Sartori, S., Trevisani, L., Nielsen, I., Tassinari, D., Ceccotti, P. and Abbasciano, V. (2003), Longevity of silicone and polyurethane catheters in long-term enteral feeding via percutaneous endoscopic gastrostomy. *Alimentary Pharmacology & Therapeutics*, 17: 853-856. <https://doi.org/10.1046/j.1365-2036.2003.01538.x>

## **POLYURETHANE BALLOON**

### **COMPLETE KIT FOR REPLACEMENT**

#### **RELIABLY REASSURING - GIVING YOU BACK TIME TO CARE**

- 1 **Distinctive white body** to easily identify as polyurethane and PURFIT™.
- 2 **Unique doughnut shaped balloon** with larger surface area internally to help reduce risk of leakage and stoma irritation.
- 3 **Inverted shaft** helps to prevent gastric wall irritation, erosion or trauma inflation and deflation.
- 4 **Tilted Inflation Valve** – to provide comfortable and easy access away from skin for inflation and deflation.
- 5 **Radiopaque stripe** along the entire length of shaft
- 6 **Balloon volume** printed on inflation valve hub.
- 7 **Purple** enteral access port with **size engraved**.
- 8 Simple and secure **press in tethered cap**.



## **SIMPLE AND EASY USER EXPERIENCE**

### **RELIABLY REASSURING**

Designed with human factors at the forefront to provide confidence and superior user, patient experience.



- ✓ **Simple press in design** – no need to marry up difficult to see markings or half twists, simply press in and hear an audible click for secure connection
- ✓ **Audible click** for extra reassurance
- ✓ **360° rotation** – use and access without restrictions
- ✓ Press **side wings** and twist for simple and easy disconnection
- ✓ Available in both **single port and dual Y-port** designs with ENFit® connectors
- ✓ **Right angle connectors** for patient comfort




## **GASTROSTOMY ACCESSORIES**

### **GASTROSTOMY UNIVERSAL CONNECTORS**

Gastrostomy Tubes are designed to be used for long term feeding and whether your patient has a low profile button device or a PEG device, it is important to maintain the functionality of the feeding tube to prevent unnecessary interruptions in Enteral Nutrition therapy or medication delivery.

To assist with the continuous use and maintenance of the device, **Hoist Medical** provide a range of universal accessories to be used with any button or PEG feeding tube. All connectors and clamps coloured purple to identify Enteral Use Only.



Universal extension sets for Low Profile Buttons. Available in single and dual Y port designs.



3 Way Stop Cock designed to connect additional ENFit® lines or syringes.



Enteral Cleaning Brush - Helps reduce accumulation and bacterial colonisation in connectors

# PRODUCT CODES



## NASO/ORAGASTRIC FEEDING TUBES INFANT & NEONATE FEEDING TUBE POLYURETHANE, RADIOPAQUE STRIPE, CLOSED DISTAL TIP

Product code	Size/description	Unit of issue
NG04-50	4Fr 50 cm	30 per box
NG04-80	4Fr 80 cm	30 per box
NG04-120	4Fr 120 cm	30 per box
NG05-50	5Fr 50 cm	30 per box
NG05-80	5Fr 80 cm	30 per box
NG05-120	5Fr 120 cm	30 per box
NG06-50	6Fr 50 cm	30 per box
NG06-80	6Fr 80 cm	30 per box
NG06-120	6Fr 120 cm	30 per box

## ENTERAL ENFit® SYRINGE PUMP

Product code	Size/description	Unit of issue
PU01	Enteral Syringe Pump	1

## EXTENSION SETS

Product code	Size/description	Unit of issue
NGE-75	75 cm	40 per box
NGE-150	150 cm	25 per box

## SYRINGES STARFIT™ ENTERAL ENFit® SYRINGES SINGLE USE

Product code	Size/description	Unit of issue
SE01LD	1 ml low dose	100 per box
SE025LD	2,5 ml low dose	100 per box
SE05LD	5 ml low dose	100 per box
SE10	10 ml	100 per box
SE20	20 ml	80 per box
SE60	60 ml	60 per box

## STARFIT™ ENTERAL ENFit® SYRINGES REUSABLE

Product code	Size/description	Unit of issue
SER01LD	1 ml low dose	100 per box
SER025LD	2,5 ml low dose	100 per box
SER05LD	5 ml low dose	100 per box
SER10	10 ml	100 per box
SER20	20 ml	80 per box
SER60	60 ml	60 per box

## ENFit® CAPS

Product code	Size/description	Unit of issue
MC01	Male Cap	100 per box
MC50	Male Cap X 50	1000 per box
MFCAP	Male Female Cap	100 per box

## MILK STRAWS

Product code	Size/description	Unit of issue
MS01	Milk Straw	50 per box
MS01WB	Milk Straw Wide Bore	50 per box

## MEDICINE DELIVERY

Product code	Size/description	Unit of issue
MFS01	Medicine Dose Straw	100 per box
MF01	Medicine Dose Aid	50 per box

## MEDICINE ACCESSORIES BOTTLE ADAPTORS

Product code	Size/description	Unit of issue
BA-L01	20-21,5 mm	100 per box
BA-M01	18,5-20,5 mm	100 per box
BA-S01	18-19,5 mm	100 per box
BA-XL01	25,5-27 mm	100 per box
BA-XS01	15-16,5 mm	100 per box
UBA01	Universal 14,5-27,5 mm	30 per box

# PRODUCT CODES



## SAMPLING STRAW

Product code	Size/description	Unit of issue
SS01	Sampling Straw	100 per box

## NEEDLE STRAW

Product code	Size/description	Unit of issue
NS01	Needle Straw	100 per box

## MALE MALE TRANSFER CONNECTOR

Product code	Size/description	Unit of issue
TR01	Male Male Transfer Connector	100 per box

## RECTAL STRAW

Product code	Size/description	Unit of issue
RK-75	Rectal Straw 75 cm	100 per box
RK-150	Rectal Straw 150 cm	50 per box

## 3 WAY STOPCOCK ENFit®

Product code	Size/description	Unit of issue
SK01	3 Way Stopcock	100 per box

## ENTERAL CLEANING BRUSH

Product code	Size/description	Unit of issue
BR01	Enteral Cleaning Brush	10 x 10

## ENTERAL DRAINAGE BAG SYSTEM DRAINAGE BAG, EXTENSION TUBING AND HANGER

Product code	Size/description	Unit of issue
DBS-100	100 ml	30 per box
DBS-250	250 ml	30 per box

## GASTROSTOMY PURFIT™ LOW PROFILE GASTROSTOMY BUTTON POLYURETHANE

Product code	Size/description	Unit of issue
PGB12-0.6	12Fr 0.6 cm	1 per box
PGB12-1.0	12Fr 1.0 cm	1 per box
PGB12-1.3	12Fr 1.3 cm	1 per box
PGB12-1.6	12Fr 1.6 cm	1 per box
PGB12-2.0	12Fr 2.0 cm	1 per box
PGB12-2.3	12Fr 2.3 cm	1 per box
PGB12-2.6	12Fr 2.6 cm	1 per box
PGB12-3.0	12Fr 3.0 cm	1 per box

## CLICK & CONNECT™ EXTENSION SETS RIGHT ANGLE, SINGLE AND DUAL Y-PORTS, STERILE

Product code	Size/description	Unit of issue
PGE-30	30 cm	6 bags of 5
PGE-60	60 cm	6 bags of 5
PGEY-30	30 cm	6 bags of 5
PGEY-60	60 cm	6 bags of 5

## UNIVERSAL BUTTON EXTENSION SETS SINGLE & DUAL Y-PORT ENFit® CONNECTORS

Product code	Size/description	Unit of issue
UGE-30	30 cm	6 bags of 5
UGE-60	60 cm	6 bags of 5
UGEY-30	30 cm	6 bags of 5
UGEY-60	60 cm	6 bags of 5

## ORDERING INFORMATION

---

### HOIST MEDICAL

HEAD OFFICE

Lot 1 - Village d'entreprise de St.Henri

Rue Anne Gacon

13016 Marseille

### ORDERS AND PRODUCT INFORMATION

Sales and General Enquiries Telephone; +33(0) 4 91 76 73 29

Department Telephone: +33(0) 6 77 74 24 74

### Orders;

[commercial@hoist-medical.com](mailto:commercial@hoist-medical.com)

### UNITED KINGDOM

[commercialuk@hoist-medical.com](mailto:commercialuk@hoist-medical.com)

[hoist-medical.com](http://hoist-medical.com)



## SUSTAINABILITY PROMISE

---

In a context where ecology has become a major concern, Hoist Medical is no exception to this trend; it is no longer just about corporate social responsibility, but a paradigm shift. Environmental protection is now at the heart of our concerns both in our daily life and in the professional context.

Global warming, pollution, deterioration of ecosystems and biodiversity, reduction of natural resources, accumulation of waste are therefore determinants of health. The challenges are numerous. The implementation of concrete actions in terms of environmental protection is part of a collective public-private awareness of the fragility of our planet's resources and of the need to develop effective policies to encourage the necessary transformations within organisations.

In addition to respecting the legislative framework, Hoist Medical implements voluntary environmental initiatives. Having become aware that the entire economy depends directly or indirectly on the good management of natural resources and the environment, as well as on biodiversity and the proper functioning of ecosystems, these approaches require prior recognition of the effects of its activity on the environment.


In this context:

- ✓ Hoist Medical is committed to economical management of resources, such as raw materials or water.
- ✓ Hoist Medical is committed to implementing steps to improve energy efficiency and reduce greenhouse gas emissions.
- ✓ Hoist Medical is committed to having an action plan in the areas of improving energy efficiency or reducing greenhouse gas emissions and the development of renewable energies.

Hoist Medical undertakes to implement actions to:

- ✓ Limitation of waste production
- ✓ Recycling, reduction and recovery of waste
- ✓ Recycling obsolete production equipment
- ✓ Integration of recycled products into these production processes
- ✓ Closed circuit water treatment or reuse
- ✓ Make eco-responsible purchases
- ✓ Support eco-gestures
- ✓ Strengthen the company culture in the area of environmental preservation

Bernard Daurelle  
President





**Hoist**  
medical

EVERY AGE, EVERY STAGE NUTRITION MATTERS



**Hoist**  
medical

**HEAD OFFICE:**

Lot 1 - Village d'entreprise de St-Henri  
Rue Anne Gacon  
13016 Marseille

Sales and General Enquiries  
+33(0) 4 91 76 73 29  
[commercial@hoist-medical.com](mailto:commercial@hoist-medical.com)

[hoist-medical.com](http://hoist-medical.com)

**UNITED KINGDOM:**

[commercialuk@hoist-medical.com](mailto:commercialuk@hoist-medical.com)

[hoist-medical.com](http://hoist-medical.com)

If you would like to know how to buy  
**STARFIT™** and **PURFIT™** products  
in your country or are interested in  
partnering with us contact  
[commercial@hoist-medical.com](mailto:commercial@hoist-medical.com)